



Health Industry Business
Communications Council

2525 E. Arizona Biltmore Cir.
Suite 127
Phoenix, AZ 85016
602-381-1091
FAX 602-381-1093
www.hibcc.org



HIBCC UDI Label Examples

DataMatrix

CompuHyper GlobalMed®

Ultra Implantable™

Fictitious Medical Device

2.25 mm x 8 mm

CAT 123ABC

LOT 1234AB



USE BY:
2019-05-15



MANUFACTURED ON:
2015-10-01



SN 5678EDFG

QTY: 1 EA



SINGLE USE



**DO NOT USE
IF PACKAGE
IS DAMAGED**



**40°C UPPER
LIMIT OF
TEMPERATURE**



KEEP DRY



Manufacturer

CompuHyper GlobalMed®
123 Technology Dr
Somewhere, XX 00000

800.555.1234 (USA)
555.555.1234 (All Others)
www.chgm.com

MedDevFront UK
Somewhereshire
XX12 3XX UK



www.mdfco.uk

+X999123ABC0/\$\$31905151234AB/S5678EDFG/16D20151001J

Linear Concatenated

CompuHyper GlobalMed®

Ultra Implantable™

Fictitious Medical Device
2.25 mm x 8 mm

CAT 123ABC

LOT 1234AB



USE BY:
2019-05-15



MANUFACTURED ON:
2015-10-01

SN 5678EDFG

QTY: 1 EA



SINGLE USE



DO NOT USE
IF PACKAGE
IS DAMAGED



40°C UPPER
LIMIT OF
TEMPERATURE



KEEP DRY

Manufacturer

CompuHyper GlobalMed®
123 Technology Dr
Somewhere, XX 00000

800.555.1234 (USA)
555.555.1234 (All Others)
www.chgm.com

MedDevFront UK
Somewhereshire
XX12 3XX UK

www.mdfco.uk



+X999123ABC0/\$\$31905151234AB/S5678EDFG/16D20151001J

Linear Not Concatenated

CompuHyper GlobalMed®

Ultra Implantable™

Fictitious Medical Device

2.25 mm x 8 mm

CAT 123ABC

LOT 1234AB



USE BY:
2019-05-15



MANUFACTURED ON:
2015-10-01

SN

5678EDFG

QTY: 1 EA



SINGLE USE




DO NOT USE
IF PACKAGE
IS DAMAGED



40°C UPPER
LIMIT OF
TEMPERATURE



KEEP DRY

 **Manufacturer**

CompuHyper GlobalMed®
123 Technology Dr
Somewhere, XX 00000

800.555.1234 (USA)
555.555.1234 (All Others)
www.chgm.com

MedDevFront UK
Somewhere
XX12 3XX UK

www.mdfco.uk



X999123ABC0B



\$S3150615:234AB/55878EDFG/16D20151001BK